

The XVu Upgrade

ADVANCING TO A NEW ERA IN HEMODYNAMIC MONITORING

XVu Comprehensive Software Upgrade.

Now is your opportunity to upgrade your present system to the new Horizon XVu. The package provides the new software features:

- ⇒ XVu Coronary diagram dedicated for graphical representation of Interventional Procedures. Editing capabilities provide easy marking of lesions, balloon, stent, graft and collateral blood vessels.
- ⇒ XVu Pullback procedure enables valve and gradient results to be displayed in the most comprehensive mode. Thumbnails provide quick access to the previously acquired waveforms.
- ⇒ XVu Full Disclosure playback panel is easily controlled. Specific patient events can be marked for quick access during review sessions.
- ⇒ XVu display layout recommended for coronary procedures where PCI table, coronary editor and AO / LV analyses are performed. Also available is a layout with a larger secondary waveform display.
- ⇒ XVu display layout recommended for Right and Left Heart procedures where detailed site analyses, heart editor and oximetry analyses are routinely performed.
- ⇒ XVu display layout recommended for procedures involving significant documentation.
- ⇒ XVu pre-configured subgroups of notes enables easy documentation of syndromes, diagnoses, procedures, inventory, etc., thus enabling detailed reporting with minimal user effort.
- ⇒ XVu PCI table allows easy and handy documentation of percutaneous interventions for review during procedures, as a table in a final report, or simply as a printed record for post review

XVu Comprehensive Hardware Upgrade

The package provides a cost effective upgrade path for the hardware and service requirements.

- ⇒ New Sun Computer with Solaris 10 operating system.
- ⇒ CFE Signal Acquisition unit upgrade.
- ⇒ Option for EtCO₂ Monitoring (system option)
- ⇒ Back-up Mode operation (system option)
- ⇒ Intra Cardiac ECG Lead Capability (6 Channels) (system option)