



Wireless Statement

“Wireless Networking” for anywhere you need patient care

The wireless network option for VitaLogik portable monitors extends central station monitoring and automated vital sign collection beyond the bedside, because “patient care” can be anywhere—at the bedside, in the hallway, in transport or in ancillary units.

Benefits include:

Utilizes an industry-standard 802.11b/g IEEE compliant radio card.

Dual-diversity dipole antenna system ensures optimal signal reception.

Designed to operate on existing wireless infrastructures that support 802.11b/g, thereby maximizing existing wireless infrastructure while managing additional costs.

Seamless network communication, whether your VL4000 or VL6000 is connected via hardwired Ethernet cable or communicating wirelessly across the hospital’s wireless local area network (WLAN).

On-screen signal strength indicator provides a clear visual of wireless networking status.

Standard WPA security and encryption to help hospitals protect patient information.

* Not available in all geographies. Please check with your local Mennen Medical representative.