

CerebraLine

New Era of Video aEEG Monitoring



The new Cerebraline offers innovative approach to the aEEG

Practice Includes:



3 Measured and 1 Calculated EEG/aEEG Channels
Simultaneous Patient Video
Quick Event Marking & Data Export
Both Built-in and External Viewers
Electrode Placement Support

A modern user interface with live EEG signals and a bold presentation of the aEEG background

Easy Events marking

Real time and post event marking – the system will save two minutes of video along with the EEG traces and aEEG background.

A special feature allows real time export of data sections during monitoring.



Live Video & EEG/aEEG

Screen Displays simultaneously patient video & EEG/aEEG traces.

Cutting edge user interface that supports various dynamic screen settings for video and/or EEG/aEEG display.

Playback screen (Built in Viewer)

Just pick the time with your fingers and the system will present the data in three simultaneous presentations of EEG, aEEG and Video.



Export Data & External Viewer on PC
The whole case can be exported and reviewed on an external Viewer (on PC)



Electrodes Placement Assistance Screen

User Interface that guides the user to pick channel locations and confirm placements.

The system uses 5 leads for three measured channels and one calculated channel.



Many NICU's in many hospitals around the world are facing the same challenge of immediate access to neurological experts.

EEG/aEEG monitoring has already gained reputation as the best solution to address this challenge.

The CerebraLine is our latest innovative addition to a vast offering of products Mennen Medical have provided to NICU's for many years.



Israel: Mennen Medical Ltd. 6 Ha-Kishon Street, Yavne, P.O. Box 102, Rehovot 7610002, Israel. Phone: +972-8-9323333. Fax: +972-8-9328510

United States: Mennen Medical Corp. 290 Andrews Road, Feasterville-Trevose, PA 19053-3480, USA. Phone: +215 259-1020. Fax: +215 357-2010

United Kingdom: Charter Kontron Ltd. Unit 18 Avant Business Centre 21 Denbigh Road Milton Keynes MK1 1DT, UK. Phone: +44 1908 646070. Fax: +44 1908 646030