


Menntor X9 Display Options & Patient Data Review

Menntor X9

Innovating Your Bedside.
Advanced configured flat-screen monitor!



Waveform + Trend Display



Remote Bed View and Alarm Watch




Full Disclosure of ECG waveforms
10 second display




Overview of Vital Signs waveforms
synchronized to top displayed ECG



Vital signs chart with up to
one minute resolution



CerebraLogik aEEG Display



We Follow Every Beat of Your heart

Israel: Mennen Medical Ltd.
6 Ha-Kishon Street, Yavne.
P.O. Box 102, Rehovot 7610002, Israel.
Phone: +972-8-9323333
Fax: +972-8-9328510
www.mennenmedical.com

United States: Mennen Medical Corp.
290 Andrews Road, Feasterville-Trevoze, PA
19053-3480, USA.
Phone: +215-259-1020
Fax: +215-357-2010
www.mennenmedical.com

United Kingdom: Charter Kontron Ltd.
Unit 18 Avant Business Center 21
Denbigh Road Milton Keynes MK1 1DT, UK.
Phone: +44-1908-646070
Fax: +44-1908-646030
www.charter-kontron.com


The latest state-of-the-art vital sign monitor

Mennen Medical is proud to introduce the newest member of its Monitoring family, the Menntor X9:

- Provides full monitoring capabilities with a state of the art 19" flat touch screen display.
- The same easy to use interface of Mennen Medical's vital signs monitoring.
- Slim all-in-one monitor for space saving and minimal/easy install into clinical areas

Features & Benefits:

- Long-term data storage:
 - Graphic Trend
 - Numeric Charts
 - ECG & Vital Signs Waveform Full Disclosure
- LAN connectivity / Wireless LAN Option
- Interface to Nurse Central Station
- Full Connectivity with Hospital network
- Easily customizable for optimal clinical use
- User friendly, intuitive interface
- No fan for a quiet environment
- **Special features** - CRG, Alarm Watch,
- **CerebraLogik line** - Dual Channel EEG & aEEG
- **Multiple Interface Options** -Anesthesia machine, PiCCO, BIS
- **Customizing Tools** - Monitor Profile, "Quick keys", Permission Editor



Big numbers display



Standard display


Menntor X9: For any patient, For every Department

The Menntor X9 adapts and fits all hospital departments and clinical needs by a wide range of vital signs monitoring, connectivity options to hospital network and vast interactive storage.

It supports Mennen Medical's Ensemble nurses station and the aEEG CerebraLogik.

It has 3 to 6 hours of battery for continuous patient monitoring and saved vital signs, for the best patient care available.

State-of-the-art 19" flat display



Optional wall mount