


Vathin Single-use Flexible Ureteroscope



Digital Video Monitor



Digital video processor



ID:
1.2mm/3.6Fr
OD:
2.8mm/8.4Fr

Specifications

Model	Working length(mm)	Work Channel ID (mm/Fr)	Outer diameter of main hose OD (mm/Fr)	Bending Range	Field of view	Bending lock	Suction function
US-S170	700	1.2 / 3.6Fr	2.8 / 8.4Fr	U285° / D285°	120°±10%	Yes	Yes