

ArcAngel

Uncuffed ET, Gastric Suction & Feeding Tube Holder

Invented by a Nurse Practitioner, ArcAngel is a secure, adjustable, multi-functional neonatal and pediatric endotracheal tube holder that eliminates tape and increases ETT/NGT placement accuracy. Silicone cheek pads are gentle, easy to clean and replaceable.




One Clinician, One Hand – Seconds for Adjustments

Key Features & Benefits.

- ✓ Securely holds both endotracheal and gastric suction/feeding tube
- ✓ Adjust ET and NG or gastric suction tubes independently of each other
- ✓ ETT/NG/gastric suction tubes can be inserted and secured both orally and nasally
- ✓ 0.25cm hash marks on the holder make adjustments safe and precise
- ✓ MRI safe – metal free
- ✓ Non allergenic – latex free
- ✓ Eliminates tape on baby's face for better visibility
- ✓ Decreased risk of skin damage due to tape irritation and repeated cleaning/reapplication
- ✓ Cheek pads can be wiped clean and are available in 2 sizes for a custom fit
- ✓ Cheek pads are designed to eliminate the need to trim
- ✓ The arc is easily removed in emergency situations without need to remove cheek pads

ArcAngel Twist Lock Technology.



Tube Holders: A. Feeding & Gastric Suction B. Endotracheal C. Feeding & Gastric Suction

Product Specifications

Size	SKU	ETT	ETT O.D. Fitment Range	Feeding & Gastric Suction Tube Sizes	
				A	C
0	WNP-AA-00	2.5mm*	3.3 - 3.7 mm	5F	6F**
1	WNP-AA-01	3.0mm	4.0 - 4.5 mm	6.5F	8F**
1x***	WNP-AA-1X	3.0mm	4.0 - 4.5 mm	6.5F	10F
2	WNP-AA-02	3.5mm	4.5 - 5.0 mm	6.5F	10F
3	WNP-AA-03	4.0mm	5.0 - 5.5 mm	6.5F	10F

* 5F inline suction catheter is recommended.

*** Same size arc as Size 2.

** Will secure 6F & 8F EDI Catheters.

Two sizes of Cheek Pads are available: Small WNP-CP-00 and Large WNP-CP-01



Safer Survival Tools for Little Warriors.



Scan here for more information

Jessica Harnish MSN, NNP-BC
President, Founder & Inventor

📞 (586) 212-3418

✉️ jessica@warriornp.com

2211 South Telegraph Road, Suite #7230
Bloomfield Hills, MI 48302